

Release Notes

Samsung Health Sensor SDK

Release Date

September 22, 2024

Release Version

v1.3.0

Target Device

The Samsung Health Sensor SDK's APIs are supported on the following watch devices that run on Wear OS Powered by Samsung.

- Galaxy Watch4 series and later models

SDK Content

Content	Content description
Documents	API Reference Programming Guide
Library	Samsung Health Sensor SDK's library. Import it in your app project. - <code>samsung-health-sensor-api.aar</code>
Sample apps	Measure blood oxygen Measure blood oxygen and heart rate Transfer measured heart rate to a connected phone Measure skin temperature

Features

The tracking APIs of the Samsung Health Sensor SDK enable partner application to register an event for tracking health sensor data and to retrieve that data while tracking on a Galaxy Watch.

Capabilities

The Samsung Health Sensor SDK provides the available tracker types on the watch.

Measuring Watch's Health Sensor Data

A partner application using the SDK can measure health-sensor data of a Galaxy Watch. The following tracker types are supported:

[Continuous Tracker Types]

The following tracker types can be measured continuously until unsetting the tracker type's event listener. Continuous tracker types operate with low battery consumption of the Galaxy Watch.

- Accelerometer
- Heart rate including IBI (inter-beat interval)
- Photoplethysmogram (PPG) green, Infrared (IR), and red
- Skin temperature

[On-demand Tracker Types]

The following tracker types are on-demand tracker types. Use only one on-demand tracker type at a time and only when needed.

- Bioelectrical impedance analysis (BIA)
- Electrocardiogram (ECG)
- PPG green, IR, and red
- Skin temperature
- SpO2 (blood oxygen)

[Other]

Measuring the user's sweat loss amount after a running is available with the following tracker type.

- Sweat loss

Developer Mode

A developer mode is supported for testing and debugging. See the developer guide for more information.

Limitations

- The emulator is not supported.
- Data measured by the Samsung Health Sensor SDK is for fitness and wellness information only, not for the diagnosis or treatment of any medical condition.

Changes

September 22, 2024

- Rebranding the SDK name to "Samsung Health Sensor SDK".

1.3.0 - August 22, 2024

[New]

- Continuous and on-demand enum constant names of `HealthTrackerType` have been defined.
 - `ACCELEROMETER_CONTINUOUS`
 - `BIA_ON_DEMAND`
 - `ECG_ON_DEMAND`
 - `HEART_RATE_CONTINUOUS`
 - `PPG_CONTINUOUS`
 - `PPG_ON_DEMAND`
 - `SKIN_TEMPERATURE_ON_DEMAND`
 - `SPO2_ON_DEMAND`
- An application can track one or more `PpgTypes` with `PPG_CONTINUOUS` or `PPG_ON_DEMAND`. The following API has been added.
 - `HealthTrackingService.getHealthTracker(HealthTrackerType, Set<PpgType>)`
- The following PPG types have been added.
 - `PpgType.GREEN`
 - `PpgType.IR`
 - `PpgType.RED`
- `ValueKey.PpgSet` has been added. It includes PPG Green, IR, Red data.
- The following fields for BIA's raw data have been added in `ValueKey.BiaSet`.
 - `BODY_IMPEDANCE_DEGREE`
 - `BODY_IMPEDANCE_MAGNITUDE`

[Deprecated]

- The following enum constant names of `HealthTrackerType` have been deprecated.
 - `ACCELEROMETER`
: Use `ACCELEROMETER_CONTINUOUS`.
 - `BIA`
: Use `BIA_ON_DEMAND`.
 - `ECG`
: Use `ECG_ON_DEMAND`.

- HEART_RATE
: Use HEART_RATE_CONTINUOUS.
 - PPG_GREEN
: Use PPG_CONTINUOUS.
 - PPG_IR
: Use PPG_ON_DEMAND.
 - PPG_RED
: Use PPG_ON_DEMAND.
 - SKIN_TEMPERATURE
: Use SKIN_TEMPERATURE_ON_DEMAND.
 - SPO2
: Use SPO2_ON_DEMAND.
- The following PPG data set of ValueKey have been deprecated. Use ValueKey.PpgSet instead of them.
 - ValueKey.PpgGreenSet
 - ValueKey.PpgIrSet
 - ValueKey.PpgRedSet

[Removed]

- The following fields of ValueKey.EcgSet have been removed. Use alternative fields.
 - ValueKey.EcgSet.ECG
: Use ValueKey.EcgSet.ECG_MV.
 - ValueKey.EcgSet.MAX_THRESHOLD
: Use ValueKey.EcgSet.MAX_THRESHOLD_MV.
 - ValueKey.EcgSet.MIN_THRESHOLD
: Use ValueKey.EcgSet.MIN_THRESHOLD_MV.
- The following fields of ValueKey.HeartRateSet have been removed. Use alternative fields.
 - ValueKey.HeartRateSet.HEART_RATE_IBI
: Use ValueKey.HeartRateSet.IBI_LIST.
 - ValueKey.HeartRateSet.STATUS
: Use ValueKey.HeartRateSet.HEART_RATE_STATUS and ValueKey.HeartRateSet.IBI_STATUS_LIST.

1.2.0 - August 30, 2023**[New]**

- The skin temperature tracker type and data point set have been added. Skin temperature can be measured on-demand or as a batching event. For more information, see the API Reference.
 - HealthTrackerType.SKIN_TEMPERATURE
 - HealthTrackerType.SKIN_TEMPERATURE_CONTINUOUS
 - ValueKey.SkinTemperatureSet

- The following interfaces for values in millivolts have been added to the ECG data point set:
 - ECG_MV
 - MAX_THRESHOLD_MV
 - MIN_THRESHOLD_MV
- The following interfaces for heart rate status and IBI-related information have been added to the heart rate data point set:
 - ValueKey.HeartRateSet.HEART_RATE_STATUS
 - ValueKey.HeartRateSet.IBI_LIST
 - ValueKey.HeartRateSet.IBI_STATUS_LIST
- The status interface has been added to the following data point sets:
 - ValueKey.AccelerometerSet.STATUS
 - ValueKey.PpgGreenSet.STATUS
 - ValueKey.PpgIrSet.STATUS
 - ValueKey.PpgRedSet.STATUS

[Change]

- The error description of the first bit in the 1 ~ 127 value for ValueKey.SweatLossSet.STATUS has been updated.

[Deprecated]

- The following interfaces in the ECG data point set have been deprecated.
 - ValueKey.EcgSet.ECG, replaced by ValueKey.EcgSet.ECG_MV.
 - ValueKey.EcgSet.MAX_THRESHOLD, replaced by ValueKey.EcgSet.MAX_THRESHOLD_MV.
 - ValueKey.EcgSet.MIN_THRESHOLD, replaced by ValueKey.EcgSet.MIN_THRESHOLD_MV.
- The following interfaces in the heart rate data point set have been deprecated.
 - ValueKey.HeartRateSet.HEART_RATE_IBI, replaced by ValueKey.HeartRateSet.IBI_LIST.
 - ValueKey.HeartRateSet.STATUS, replaced by ValueKey.HeartRateSet.HEART_RATE_STATUS and ValueKey.HeartRateSet.IBI_STATUS_LIST.

The following values of ValueKey.SweatLossSet.STATUS have been deprecated.

STATUS value	Description
-10	A PPG signal is too weak. An accuracy is very low.
-8	A PPG signal is weak. An accuracy is low.
-3	A Galaxy wearable is detached.
-1	Sweat loss sensor's accelerometer data accuracy is low.

[Removed]

- `ValueKey.HeartRateSet.FLAG` has been removed.

1.1.0 - March 31, 2022

- Minor stability improvements have been made to the SDK library.

1.1.0 Alpha - March 17, 2022

[New]

- The `SpO2` tracker type and its data point set have been added. It is supported by Health Platform v1.3.0.
 - `HealthTrackerType.SPO2`
 - `ValueKey.Spo2Set`
- The `ValueKey.HeartRateSet.STATUS` has been added, which replaces the `ValueKey.HeartRateSet.FLAG`.

[Changes]

- The `TrackingSampleApp` application code has been updated.

[Deprecated]

- The `ValueKey.HeartRateSet.FLAG` has been deprecated.

1.0.0 - Feb 25, 2022

- The values of the `SweatLossSet.STATUS` have changed. For more information, see API Reference.

1.0.0 beta1 - Dec 31, 2021

- Minor issues related to the sweat loss feature have been fixed.

1.0.0 alpha3 - Dec 17, 2021

[New]

- The sweat loss tracker type has been added. It measures how much sweat is lost from a running. The following class, enums, and APIs have been added:
 - `DataType`
 - `ExerciseState`
 - `ExerciseType`
 - `HealthTrackerType.SWEAT_LOSS`
 - `HealthTrackingService.getHealthTracker(HealthTrackerType healthTrackerType, TrackerUserProfile userProfile, ExerciseType type)`
 - `ValueKey.SweatLossSet`

[Known Issue]

- The following issues are to be solved in February 2022 watch software update:
 - A timestamp period variation issue in received sensor data for `AccelerometerSet`, `HeartRateSet`, and `PpgGreenSet`.
 - The `HealthTracker.TrackerEventListener.onFlushCompleted()` is not received.

1.0.0 alpha2 - July 28, 2021**[Changes]**

- The partner application's SDK policy has been applied. Tracking watch sensor data with the SK is available only for the registered tracker type scope.
If `SDK_POLICY_ERROR` occurs, ask to your Samsung contact point.
- The following tracking data point sets have added more interfaces.
 - `ValueKey.BiaSet.STATUS`
: The BIA status has been defined. Check this value after the BIA measurement. If it indicates an error guide the user with a proper UI and a message.
 - `ValueKey.EcgSet.LEAD_OFF`
: The ECG's lead on/off values have been defined.
 - `ValueKey.HeartRateSet.FLAG`
: The flag for heart rate measurement has been added.
- In `HealthTracker.TrackerError`:
 - The following unused errors have been removed.

- LOW_SIGNAL
- MOVEMENT_DETECTED
- NOT_WEARING
- TIME_OUT
- A new error has been added.
 - SDK_POLICY_ERROR
- Measurement guides for each sensor data have been added. See Programing Guide.